

ONH836

Oxygen/Nitrogen/Hydrogen

This instrument now supports either Helium or Argon carrier gas. The type of carrier gas used may affect some instrument specifications, as indicated below.

Specification Sheet

| Instrument Range* (1 g sample) | | Precision† | |
|--|--|---------------------------|--|
| Oxygen: | 0.00005 mg** to 50 mg | Oxygen: | 0.000025 mg or 0.3 % RSD, whichever is greater |
| Nitrogen, He Carrier Gas: | 0.00005 mg** to 30 mg | Nitrogen, HE Carrier Gas: | 0.000025 mg or 0.3 % RSD, whichever is greater |
| Nitrogen, Ar Carrier Gas: | 0.0002 mg** to 30 mg | Nitrogen, Ar Carrier | 0.0001 mg or 0.3 % RSD, whichever is greater |
| Hydrogen: | 0.0001 mg** to 2.5 mg | | 0.00005 mg or 2 % RSD, whichever is greater |
| Analysis Time§ | | | |
| Oxygen, He Carrier Gas: | 85 s | Oxygen, Ar Carrier Gas: | 95 s |
| Nitrogen, He Carrier Gas: | 100 s | Nitrogen, Ar Carrier Gas: | 130 s |
| Hydrogen, He Carrier Gas: | 90 s | | 100 s |
| Cycle Time§ (including outgas, purge, analysis delay, and analysis time) | | | |
| He Carrier Gas: | 180 s | Ar Carrier Gas: | 210 s |
| Calibration | Standards (single or muöti-point); manual, gas dose | | |
| Sample Size | 1 g (nominal) | | |
| Detection Method | Non-Dispersive Infrared Absorption; Thermal Conductivity | | |
| Chemical Reagents | Anhydrous Magnesium Perchlorate (MgClO4) Copper Oxide, Rare Earth Copper Oxide,Copper Turnings | | Oxygen/Moisture Indicating Tube Sodium Hydroxide on an Inert Base |
| Gas Requirements | | | |
| He carrier Gas: | Helium (99.99 % pure), 22 psi (1.5 bar) ± 5 % | | |
| Ar Carrier Gas: | Argon (99.999 % pure), 22 psi (1.5 bar) ± 5 % | | |
| Pneumatic: | Compressed Air, 40 psi (2.8 bar) ± 10 %, source must be oil and water free | | |
| Gases Optional | | | |
| Gas Dose: | Carbon Dioxide, 99.99 % pure, 20 psi (1.4 bar) ±10 % | | |
| Gas Dose: | Nitrogen, 99.99 % pure, 20 psi (1.4 bar) ± 10 % | | |
| Gas Flow Rates | Carrier: 490 cm³/min | | Pneumatic: 280 cm³/analysis |
| Furnace | Impulse furnace with current and power control 7500 W maximum, liquid cooled | | |
| Coolant | 3.2 L LECO Coolant | | |
| Operation Conditions | Temperature: 15 °C to 35 °C (59 °F to 95 °F) | | Rel. Humidity: 20 % to 80 %, non-condensing |
| Sound Pressure Level | Operating Temp: 15 °C to 35 °C (59 °F to 95 °F) | | Rel. Humidity: 20 % to 80 % (non-condensing) |
| Dimensions†† | 36 in H x 28 in W x 34 in D (91 cm x 71 cm x 86 cm) with touch-screen monitor | | |
| Electrical Power | 230 V~ (+ 10/-15 %; at max load); 50 A, 50/60 Hz, Single Phase; 12500 Btu/h‡ | | |
| Weight (approx.) | Analyzer: 410 lb (186 kg) without touch-screen monitor | | |

Part Numbers

| | |
|--------------|--|
| ONH836-XXXXC | Oxygen/Nitrogen Determinator with software and external PC |
|--------------|--|

Options

NOTE: Multiple configurations of options are available. Please contact your local LECO Sales Engineer for more details.

- Unit with PC, touch-screen monitor package (C)
- Optional mounted touch-screen monitor package (M)
- Optional automatic cleaner package (H)



- * Use the following formula to calculate element concentration:
% Element concentration = [(absolute element mass in mg)/(sample mass in mg)]*100.
- ** Lower range is calculated as 2σ instrument blank deviation. Method range may differ due to factors such as sample type and method parameters.
- † Calculated as 1σ instrument blank deviation. Method precision may differ due to sample inhomogeneity or other external factors.
- § All times listed are nominal, actual times may vary based on method settings and application.
- †† Allow for a 6 in (15 cm) minimum access area around all sides.
- ‡ Average output based on nominal operating parameters.

Theory of Operation

The ONH836 Oxygen/Nitrogen/Hydrogen system is designed for simultaneous wide-range measurement of oxygen, nitrogen, and hydrogen content of steel, refractory metals, and other inorganic materials. The instrument features custom software designed specifically for touch operation.

A pre-weighed sample is placed in a graphite crucible which is heated in an impulse furnace to release analyte gases.


Oxygen present in the sample reacts with the graphite crucible to form CO and CO₂. An inert gas carrier, typically helium, sweeps the liberated analyte gases out of the furnace, through a Mass Flow Controller, and through a series of detectors. CO and CO₂ are detected using non-dispersive infrared (NDIR) cells. The gas then flows through a heated reagent where the CO is oxidized to form CO, and H is oxidized to form H₂O. The gas continues through another set of NDIR cells where H₂O and CO₂ are detected. The CO and H₂O are then scrubbed out of the carrier gas stream, leaving the final analyte, nitrogen, as the only impurity. A patented Dynamic Flow Compensation (DFC) system is used to add carrier gas as a makeup for the gas lost during the scrubbing process. A Thermal Conductivity (TC) detector is used to detect nitrogen.

The detection system is comprised of both NDIR and TC detectors. NDIR cells are based on the principle that analyte gas molecules absorb infrared (IR) energy at unique wavelengths within the IR spectrum. Incident IR energy at these wavelengths is absorbed as the gases pass through the IR absorption cells. The complete set of CO and CO₂ NDIR cells is required to give the most accurate oxygen results for a wide range of sample types and concentrations. TC detection takes advantage of the difference in thermal conductivity between carrier and analyte gases. Resistive TC filaments are placed in a flowing stream of carrier gas and heated by a bridge circuit. As analyte gas is introduced into the carrier stream, the rate at which heat transfers from the filaments will change producing a measurable deflection in the bridge circuit.

The concentration of an unknown sample is determined relative to calibration standards. To reduce interferences from instrument drift, reference measurements of pure carrier gas are made prior to each analysis.

Flow Diagram Shown for ONH836 Configuration

Configuration-specific flow diagrams available on request



Specifications and part numbers may change.
Consult LECO for latest information.

LECO Corporation
3000 Lakeview Avenue | St. Joseph, MI 49085
Phone: 269-985-5496
info@leco.com | www.leco.com

LECO Europe
eu.leco.com

LECO
EMPOWERING RESULTS